

Roche KAPA EvoPlus V2 and KAPA EvoPlus Boost



Overview

Roche KAPA EvoPlus Boost Kits and KAPA EvoPlus V2 Kits are ideally suited for low- and high-throughput Next-Generation Sequencing (NGS) library construction workflows requiring enzymatic fragmentation. These kits are designed for library construction from a wide range of sample types and inputs (0.1 ng to 500 ng) and are automation friendly.

Full walkway KAPA EvoPlus V2 (or KAPA EvoPlus Boost) protocols have been developed for the SPT Labtech firefly+. This application note provides supporting information for the firefly+ protocols listed below. These protocols are based on the Roche document - "Instructions For Use of KAPA DNA Library Prep Evolved Workflows with enzymatic fragmentation. For use with the KAPA EvoPlus V2 and KAPA EvoPlus Boost kits. Featuring KAPA EvoT4 DNA Ligase and KAPA EvoAmp ReadyMix. (October 2025, Version 5.0)" and are available to download from the firefly community.

Here, we demonstrate how these workflows have been automated on firefly+. We provide details on the steps performed in each protocol, together with data to demonstrate the performance of these protocols in generating high quality libraries for sequencing.

Protocol highlights

- Full walkway library preparation protocol for KAPA EvoPlus V2 (or KAPA EvoPlus Boost)
- Full walkway library preparation protocol for KAPA EvoPlus V2 (PCR-Free)
- On-deck thermocycler (ODTC) on firefly+ allows all incubation and thermocycling steps to take place on firefly+
- Flexibility to process 1 to 12 sample columns per run and to specify the starting column for the reagent plates (ReadyMix plates and UDI Adapters plates) - enabling multiple low-throughput runs using the same reagent plates and reducing waste.
- Additional protocols provide a flexible framework to enable multiple workflow options to be automated, to include:
 - Low DNA input (0.1ng) workflows
 - An optional double-sided size selection can be run either at the post ligation or post amplification stage
 - Library Amplification can be run with either a single or double post PCR purification

firefly+ protocols

Protocol name	firefly+ run time
KAPA EvoPlus V2 (or KAPA EvoPlus Boost)	2 hours 54 minutes (8 samples with 5 cycles PCR) to 3 hours 20 minutes (96 samples with 5 cycles PCR)
KAPA EvoPlus V2 (PCR-Free)	1 hour 56 minutes (8 samples)
KAPA EvoPlus V2 (Appendix A. Double-sided Size Selection)	42 minutes (8 samples)
KAPA EvoPlus V2 (or KAPA EvoPlus Boost) - Chapter 4. Steps 1-3a only	56 minutes (8 samples with 5 cycles PCR)
KAPA EvoPlus V2 (or KAPA EvoPlus Boost) - Chapter 4. Steps 1-3b only	1 hour 38 minutes (8 samples with 13 cycles of PCR)

Table 1. firefly+ protocol names and run times.

firefly+ protocol to run	
Workflow	KAPA EvoPlus V2 (PCR-Free)
KAPA EvoPlus V2 (PCR-Free)	
KAPA EvoPlus V2 (PCR-Free) with post ligation size selection	

Table 2. Summary of the KAPA EvoPlus V2 PCR-Free workflows that can be run and the firefly+ protocols required to run them. The KAPA EvoPlus PCR-Free workflow can be run as a standalone walkway protocol or with an additional post ligation size selection.

firefly+ protocol to run					
Workflow	KAPA EvoPlus V2 (PCR-Free)	KAPA EvoPlus V2 (or KAPA EvoPlus Boost)	KAPA EvoPlus V2 (Appendix A. Double-sided Size Selection)	KAPA EvoPlus V2 (and EvoPlus Boost) - Chapter 4. Steps 1-3a only	KAPA EvoPlus V2 (and EvoPlus Boost) - Chapter 4. Steps 1-3b only
KAPA EvoPlus V2					
KAPA EvoPlus V2 with post ligation size selection					
KAPA EvoPlus V2 with post amplification size selection					
Low input (0.1ng) KAPA EvoPlus V2					

Table 3. Summary of the KAPA EvoPlus V2 workflows that can be run and the firefly+ protocols required to run them. The KAPA EvoPlus V2 workflow uses the KAPA HiFi HotStart ReadyMix and associated PCR cycling conditions. It can be run standalone or with a post ligation or post amplification size selection or with low input (0.1ng) by running the relevant combination of protocols.

firefly+ protocol to run					
Workflow	KAPA EvoPlus V2 (PCR-Free)	KAPA EvoPlus V2 (and EvoPlus Boost)	KAPA EvoPlus V2 (Appendix A. Double-sided Size Selection)	KAPA EvoPlus V2 (and EvoPlus Boost) - Chapter 4. Steps 1-3a only	KAPA EvoPlus V2 (and EvoPlus Boost) - Chapter 4. Steps 1-3b only
KAPA EvoPlus Boost					
KAPA EvoPlus Boost with post ligation size selection					
KAPA EvoPlus Boost with post amplification size selection					
Low input (0.1ng) KAPA EvoPlus Boost					

Table 4. Summary of the KAPA EvoPlus Boost workflows that can be run and the firefly+ protocols required to run them. The KAPA EvoPlus Boost workflow uses the KAPA EvoAmp ReadyMix and associated PCR cycling conditions. It can be run standalone or with a post ligation or post amplification size selection or with low input (0.1ng) by running the relevant combination of protocols.

Protocol overview

SPT Labtech firefly+ protocols have been developed to run the KAPA EvoPlus V2 PCR-Free, KAPA EvoPlus V2 and KAPA EvoPlus Boost workflows - using the 96-well plate format kits and the KAPA Unique Dual-Indexed Adapter Kit. The firefly+ protocols run in 96 well plate format and make use of the on-deck thermocycler (ODTC) on firefly+ to maximise walkway time. Variables allow a user to run from 1 to 12 columns of library preparation. Either the KAPA EvoPlus V2 or KAPA EvoPlus Boost workflow can be run using the same firefly+ protocols.

The automated KAPA EvoPlus V2 workflow should be used with the KAPA EvoPlus V2 Kit (96-well plate). This kit contains a KAPA HiFi HotStart ReadyMix (2X) plate which needs to be used with the KAPA HiFi HotStart ReadyMix cycling program conditions - by selecting the relevant "KAPA HiFi HotStart ReadyMix cycling program" PCR ODTC protocol when starting the run.

The automated KAPA EvoPlus Boost workflow should be used with the KAPA EvoPlus Boost Kit (96-well plate). This kit contains a KAPA EvoAmp ReadyMix plate which needs to be used with the KAPA EvoAmp ReadyMix cycling program conditions - by selecting the relevant "KAPA EvoAmp ReadyMix cycling program" PCR ODTC protocol when starting the run.

The KAPA EvoPlus V2 and KAPA EvoPlus Boost workflows on firefly+ cover the following sections of the Instructions For Use document:

Chapter 3. Prepare the Sample Library:

- Step 1. Fragmentation and A-tailing
- Step 2. Adapter Ligation
- Step 3. Purify the Sample Library using KAPA HyperPure Beads

Chapter 4. Amplify the Sample Library:

- Step 1. Prepare the Library Amplification Reaction
- Step 2. Perform the Library Amplification
- Step 3. Purify the Amplified Sample Library using KAPA HyperPure Beads

Appendix A. Double-sided Size Selection

KAPA EvoPlus V2 (PCR-Free) firefly+ protocol

This full walkway protocol runs the KAPA EvoPlus V2 (PCR-Free) workflow on firefly+ following Chapter 3. Prepare the Sample Library of the Instructions for use document.

The protocol adds FragTail ReadyMix to the input DNA and tip mixes these reagents. The input DNA plate is then placed on the On Deck Thermocycler (ODTC) to run the Fragmentation and A-tailing. On completion of the Fragmentation and A-tailing, the protocol adds KAPA UDI Adapters and Ligation ReadyMix to the FragTail product and tip mixes these reagents. The plate is then incubated at 20°C for 15 minutes on the ODTC. Following the Adapter Ligation the protocol performs a 0.6X bead purification and then transfers the cleaned-up adapter ligated libraries to a fresh plate.




Figure 1. firefly+ starting deck layout for a 96-sample KAPA EvoPlus V2 (PCR-Free) library prep run. Deck layout - Basement: Reservoir R1 - Elution Buffer, R2+R3+R5+R6 - 80% Ethanol, R4 - KAPA HyperPure Beads, Plate Thermal Module - 96 sample Thermo Adapter Block and Input DNA Plate. Lower Deck: L1+L4+L5+L6+L7+L8 - 100µL 96 format ATL tips on ATL 35-125µL Tip Stands, L3 - Waste Plate. Upper Deck: U2 - FragTail ReadyMix Plate, U3 - Ligation ReadyMix Plate, U4 - KAPA UDI Adapter Plate, U7 - KAPA HyperPure Beads Plate, U8 - Alpaqua Magnum FLX 96 sample Magnetic Block.



Figure 2. firefly+ starting layout for a 96-sample KAPA EvoPlus V2 (PCR-Free) library prep run. M1-M3 - 100µL 96 format ATL tips. S1 - Cleaned up Adapter Ligation Library Plate. S7, S8 - Plate Lids.

Once the thermocycling program is complete the plate is returned to the firefly deck. The protocol then performs a 1.0X bead purification and transfers the cleaned up amplified libraries to a fresh plate.




Figure 3. firefly+ starting deck layout for a 96-sample KAPA EvoPlus V2 (and EvoPlus Boost) library prep run. Deck layout - Basement: Reservoir R1 - Elution Buffer, R2+R3+R5+R6 - 80% Ethanol, R4 - KAPA HyperPure Beads, Plate Thermal Module - 96 sample Thermo Adapter Block and Input DNA Plate. Lower Deck: L1+L4+L5+L6+L7+L8 - 100µL 96 format ATL tips on ATL 35-125µL Tip Stands, L3 - Waste Plate. Upper Deck: U2 - FragTail ReadyMix Plate, U3 - Ligation ReadyMix Plate, U4 - KAPA UDI Adapter Plate, U5 - Cleaned Up Adapter Ligation Library Plate, U7 - KAPA HyperPure Beads Plate, U8 - Alpaqua Magnum FLX 96 sample Magnetic Block.



Figure 4. firefly+ starting layout for a 96-sample KAPA EvoPlus V2 (and EvoPlus Boost) library prep run. M1-M11 - 100µL 96 format ATL tips. S1 - KAPA Library Amplification Primer Mix Plate, S2 - KAPA EvoAmp ReadyMix OR KAPA HiFi HotStart ReadyMix (2X) Plate, S3 - Cleaned up Amplified Library Plate, S7-S9 - Plate Lids.




Figure 5. firefly+ starting deck layout for a 96-sample KAPA EvoPlus V2 (Appendix A. Double-sided Size Selection) run. Deck layout - Basement: Reservoir R1 - Elution Buffer, R2+R3+R5+R6 - 80% Ethanol, R4 - KAPA HyperPure Beads. Lower Deck: L1+L4+L6+L7+L8 - 100µL 96 format ATL tips on ATL 35-125µL Tip Stands, L5 - Empty Waste tip Box. Upper Deck: U1 - Waste Plate, U2 - Intermediate Plate, U3 - Output Plate, U6 - Input Plate, U8 - Alpaqua Magnum FLX 96 sample Magnetic Block.



Figure 6. firefly+ starting layout for a 96-sample KAPA EvoPlus V2 (Appendix A. Double-sided Size Selection) run. M1 - Empty Waste tip Box, M2-M3 - 100µL 96 format ATL tips.

KAPA EvoPlus V2 (or KAPA EvoPlus Boost) firefly+ protocol

This full walkway protocol runs the KAPA EvoPlus V2 (or KAPA EvoPlus Boost) workflow on firefly+ following Chapter 3. Prepare the Sample Library of the Instructions for use document. This firefly+ protocol runs the same steps as the PCR-Free protocol and then carries on to the amplification and purification steps in Chapter 4 Steps 1-3a Purification of amplified Sample Library procedure.

This protocol conducts the same steps as described for the "KAPA EvoPlus V2 (PCR-Free) firefly+ protocol" above; then, following the Adapter Ligation Purification, this protocol adds Library Amplification Primer Mix and KAPA HiFi HotStart ReadyMix (2X) OR KAPA EvoAmp ReadyMix to the Adapter-ligated libraries and tip mixes. The ReadyMix required depends on the workflow being run. The relevant user defined KAPA Library Amplification program is then run on the ODTC. The KAPA EvoPlus V2 workflows use the KAPA HiFi HotStart ReadyMix (2X) plate and one set of cycling conditions. The KAPA EvoPlus Boost workflows use the KAPA EvoAmp plate and an alternative set of cycling conditions.

Once the thermocycling program is complete the plate is returned to the firefly deck. The protocol then performs a 1.0X bead purification and transfers the cleaned up amplified libraries to a fresh plate.

KAPA EvoPlus V2 (Appendix A. Double-sided Size Selection) firefly+ protocol

This optional Size Selection protocol can be run after the firefly+ KAPA EvoPlus V2 (PCR-Free) protocol - to perform post-ligation sample purification, or after the firefly+ KAPA EvoPlus V2 (and EvoPlus Boost) protocol - to perform a post-amplification sample purification. The input volume for the size selection is 50µL and output volume is 20µL. This protocol performs a double-sided Size Selection purification where the user selects the SPRI 1 and SPRI 2 ratios. The Size Selected libraries are then transferred to a fresh output plate.

Protocol performance

The KAPA EvoPlus V2 (with the KAPA EvoPlus Boost workflow) library preparation protocols were run on firefly+ and manually. Three tests were performed: Test 1 was a high-throughput run, using 96 samples at 10ng input, including 15 NTCs; Test 2 was a PCR-Free run at 500ng input (and then later a Size Selection) and Test 3 was a 0.1ng low input run. See Table 5 for details. The ODTC was used for thermocycling steps for both the firefly+ and manual runs.

The resulting libraries were analyzed to determine their concentration by qPCR and average fragment size by using a LightCycler 480 System (Roche, KAPA Library Quantification kit) and a Fragment Analyzer (Agilent, DNF-474 HS NGS Fragment kit) respectively. See Table 6 for a summary of the results. Figures 7-10 show a further breakdown of the results. No library was detected in the NTCs - see Figure 9.

Experiment	Input (ng)	firefly+ replicates	Manual replicates	Adapter Conc (µM)	PCR Cycles	Post Amplification cleanup
Test 1 high-throughput	10	96	8	6	5	Chapter 4 Step 3a (Single 1.0X cleanup)
Test 2 (PCR-Free)	500	8	4	15	0'	N/A
Test 3 low input	0.1	8	8	0.6	13	Chapter 4 Step 3b (Double 1.0X cleanup)

Table 5. Experiments run to assess protocol performance of KAPA EvoPlus Boost on firefly+. *5 cycles of PCR run for sizing purposes.

Experiment	Number samples sized	Average Size (bp)	Average Size CV%	Number of Samples quantified	Concentration (nM)	Concentration CV%
Test 1 high-throughput - firefly+	16	351	1.3	80	6.1	13.6
Test 1 manual + ODTC	8	347	2.6	8	7.15	13.4
Test 2 PCR-Free (Post Ligation) - firefly+	8	332	1.6	8	49.78	10.5
Test 2 PCR-Free (Post Ligation) - manual	4	342	3.8	4	33.73	12.0
Test 3 low input - firefly+	8	362	2.1	8	15.67	28.6
Test 3 low input - manual + ODTC	8	394	4.5	8	11.47	27.7

Table 6. Summary of sizing and concentration data for KAPA EvoPlus Boost firefly+ and manual runs.




Figure 7. Fragment Analyzer average size (bp) results for KAPA EvoPlus Boost libraries for each of the 3 tests - 10ng input with 5 cycles of PCR, 500ng input PCR-Free and 0.1ng input with 13 cycles of PCR. The comparison between libraries generated on the firefly+ and manually is shown.




Figure 8. Fragment Analyzer traces for samples run on test 1 10ng input firefly+ n=16.




Figure 9. Fragment Analyzer traces for NTCs run on test 1 10ng input firefly+ n=4.




Figure 10. KAPA Library Quantification size adjusted concentration (nM) results for KAPA EvoPlus Boost libraries for each of the 3 tests - 10ng input with 5 cycles of PCR, 500ng input PCR-Free and 0.1ng input with 13 cycles of PCR. The comparison between libraries generated on the firefly+ and manually is shown. NTCs run as part of firefly+ test 1 are also shown.

Size Selection

As part of Test 2, the 500ng input PCR-Free samples underwent 5 cycles of PCR to access sizing, a post PCR cleanup and then an additional double-sided Size Selection cleanup following the post PCR cleanup, with the first cut at 0.5X and the second cut at 0.7X. Library Size Pre and Post Size Selection was measured for the firefly+ and manual runs using Fragment Analyzer and the results are shown in Table 7 and Figures 11-12. The results show that a double-sided size selection was able to be run successfully on firefly+ with comparable results to manual size selection.

Run type	Average Size (bp) Pre Size Selection	Average Size (bp) Post Size Selection	Number of Samples
firefly+	368	429	8
Manual	362	408	3

Table 7. Summary of the Fragment Analyzer sizing data for the KAPA EvoPlus V2 (EvoAmp) libraries undergoing double sided size selection.



Figure 11. Fragment Analyzer average size (bp) for KAPA EvoPlus V2 (EvoAmp) libraries undergoing double-sided size selection on firefly+ and run manually.




Figure 12. Test 2 Size Selection on firefly+. Fragment Analyzer traces shown are for sample in position A1 pre and post size selection. Size selection was performed following a 5 cycle PCR and cleanup. 1 in 10 dilution run on the fragment analyzer.

Reagents - Kits

Description	Roche Cat No.
KAPA EvoPlus V2 Kit, plated format (96wn)	0942042001
KAPA EvoPlus V2 Kit, PCR-free, plated format (96wn)	0942043001
KAPA Library Amp Primer Mix 96-well plate (96wn)	0942047001
KAPA EvoPlus Boost Kit, plated format (96wn)	1061372001
KAPA Unique Dual-Indexed Adapter Kit, (15 µM)	0886191972
KAPA Adapter Dilution Buffer (25 mL)	0827833001
KAPA HyperPure Beads Kit (60mL)	8963851001

Summary

These results demonstrate that the KAPA EvoPlus V2 (KAPA EvoPlus Boost) and KAPA EvoPlus V2 (PCR-Free) workflows have been successfully automated on firefly+ to generate illumina-compatible sequencing libraries. Full walkway firefly+ protocols are available for KAPA EvoPlus V2 (and EvoPlus Boost) and KAPA EvoPlus V2 (PCR-Free). Additional firefly+ protocols are available to perform optional double-sided size selection and alternative versions of the library amplification and cleanup steps. Plate-based reagents and user input variables in the firefly+ protocols allow flexibility to run varying numbers of columns of samples and to start at varying columns in the reagent plates - allowing both high-throughput or multiple low-throughput run approaches from a single kit. Libraries prepared on firefly+ are in line with manually prepared libraries in terms of average fragment size and concentration. The libraries show consistent fragment sizing and concentration across a 96-well plate and show no detectable well-to-well contamination. The firefly+ size selection workflow also shows expected results consistent with manually prepared libraries.