



## firefly<sup>®</sup> technical note

# Kinnex full-length RNA Kit Library Prep

This technical note provides supporting information for automating PacBio's Kinnex full-length RNA Kit (103-238-700 REV06) Library Prep on SPT Labtech firefly liquid handler. These protocols are available to download from the firefly community. Here, we outline protocol run times, parts required and provide details on the steps performed in each protocol.

## firefly protocols

| Protocol number   | Protocol name                              | Estimated run time (minutes) |
|-------------------|--|------------------------------|
| Protocol 1 of 10  | 2.0A cDNA Synthesis                        | 10                           |
| Protocol 2 of 10  | 2.0B 1.3X SMRTbell bead cleanup            | 35                           |
| Protocol 3 of 10  | 3.0A cDNA amplification                    | 5                            |
| Protocol 4 of 10  | 3.0B 0.9X cDNA SMRTbell bead cleanup       | 35                           |
| Protocol 5 of 10  | 4.0A Kinnex PCR                            | 15                           |
| Protocol 6 of 10  | 4.0B Pooling & 1.05X SMRTbell bead cleanup | 40                           |
| Protocol 7 of 10  | 5.0A Kinnex array formation                | 10                           |
| Protocol 8 of 10  | 5.0B 1X SMRTbell bead cleanup              | 35                           |
| Protocol 9 of 10  | 6.0A Nuclease treatment                    | 5                            |
| Protocol 10 of 10 | 6.0B Final SMRTbell cleanup                | 35                           |

**Table 1.** Protocols & estimated run times used in Kinnex full-length RNA kit Library Prep on firefly.

# Input variables

| Protocol number   | Protocol name                              | Variable ID                         | Default Value |
|-------------------|--|-------------------------------------|---------------|
| Protocol 1 of 10  | 2.0A cDNA Synthesis                        | Number of Samples                   | 24            |
| Protocol 1 of 10  | 2.0A cDNA Synthesis                        | LDV Reservoir Dead Volume (μL)      | 55            |
| Protocol 2 of 10  | 2.0B 1.3X SMRTbell bead cleanup            | Sample Start Column                 | 1             |
| Protocol 2 of 10  | 2.0B 1.3X SMRTbell bead cleanup            | Number of Samples                   | 24            |
| Protocol 2 of 10  | 2.0B 1.3X SMRTbell bead cleanup            | Supernatant Removal Starting Column | 4             |
| Protocol 3 of 10  | 3.0A cDNA amplification                    | Number of Samples                   | 24            |
| Protocol 3 of 10  | 3.0A cDNA amplification                    | LDV Reservoir Dead Volume (μL)      | 55            |
| Protocol 3 of 10  | 3.0A cDNA amplification                    | Primer Plate Starting Column        | 1             |
| Protocol 4 of 10  | 3.0B 0.9X cDNA SMRTbell bead cleanup       | Sample Starting Column              | 1             |
| Protocol 4 of 10  | 3.0B 0.9X cDNA SMRTbell bead cleanup       | Number of Samples                   | 24            |
| Protocol 4 of 10  | 3.0B 0.9X cDNA SMRTbell bead cleanup       | Elution Volume (μL)                 | 24            |
| Protocol 4 of 10  | 3.0B 0.9X cDNA SMRTbell bead cleanup       | cDNA Resuspension Volume (μL)       | 55            |
| Protocol 5 of 10  | 4.0A Kinnex PCR                            | LDV Reservoir dv (μL)               | 55            |
| Protocol 5 of 10  | 4.0A Kinnex PCR                            | Number of PCR's / Kinnex Master Mix | 8             |
| Protocol 5 of 10  | 4.0A Kinnex PCR                            | Kinnex Master Mix Volume (μL)       | 22.5          |
| Protocol 5 of 10  | 4.0A Kinnex PCR                            | Kinnex Primer Starting Column       | 1             |
| Protocol 5 of 10  | 4.0A Kinnex PCR                            | Kinnex Master Mix's                 | 6             |
| Protocol 5 of 10  | 4.0A Kinnex PCR                            | Total Number of Kinnex PCRs         | 48            |
| Protocol 6 of 10  | 4.0B Pooling & 1.05X SMRTbell bead cleanup | Sample Starting Column              | 1             |
| Protocol 6 of 10  | 4.0B Pooling & 1.05X SMRTbell bead cleanup | Number of Samples                   | 16            |
| Protocol 6 of 10  | 4.0B Pooling & 1.05X SMRTbell bead cleanup | Supernatant Removal Starting Column | 3             |
| Protocol 6 of 10  | 4.0B Pooling & 1.05X SMRTbell bead cleanup | QC Aliquot Starting Column          | 3             |
| Protocol 7 of 10  | 5.0A Kinnex array formation                | Number of Samples                   | 8             |
| Protocol 7 of 10  | 5.0A Kinnex array formation                | LDV Reservoir dv (μL)               | 55            |
| Protocol 8 of 10  | 5.0B 1X SMRTbell bead cleanup              | Starting column                     | 1             |
| Protocol 8 of 10  | 5.0B 1X SMRTbell bead cleanup              | Number of samples                   | 8             |
| Protocol 8 of 10  | 5.0B 1X SMRTbell bead cleanup              | Supernatant Removal Starting Column | 2             |
| Protocol 8 of 10  | 5.0B 1X SMRTbell bead cleanup              | Elution Volume (μL)                 | 41            |
| Protocol 8 of 10  | 5.0B 1X SMRTbell bead cleanup              | QC Aliquot Starting Column          | 2             |
| Protocol 9 of 10  | 6.0A Nuclease treatment                    | Number of Samples                   | 8             |
| Protocol 9 of 10  | 6.0A Nuclease treatment                    | LDV Reservoir dv (μL)               | 55            |
| Protocol 10 of 10 | 6.0B Final SMRTbell cleanup                | Starting Column                     | 1             |
| Protocol 10 of 10 | 6.0B Final SMRTbell cleanup                | Number of Samples                   | 8             |
| Protocol 10 of 10 | 6.0B Final SMRTbell cleanup                | Supernatant Removal Starting Column | 2             |
| Protocol 10 of 10 | 6.0B Final SMRTbell cleanup                | Elution Volume (μL)                 | 26            |
| Protocol 10 of 10 | 6.0B Final SMRTbell cleanup                | QC Aliquot Starting Column          | 2             |

**Table 2.** Variables used in Kinnex full-length RNA kit Library Prep on firefly. Static variables, including those defined as algebraic expressions, are not shown.

## Reagent volumes

The reagent volumes required to run Kinnex full-length RNA kit Library Prep on SPT Labtech firefly depend on the number of samples being processed. Default required minimum volumes for these reagents, based on the number of samples shown in the **Input variables** table, are shown below and in the EXECUTE section of the firefly software.

### Protocol 1 of 10

#### 2.0A cDNA Synthesis

##### REAGENTS



Figure 1. 2.0A cDNA Synthesis minimum required reagent volumes.

### Protocol 2 of 10

#### 2.0B 1.3X SMRTbell bead cleanup

##### REAGENTS



Figure 2. 2.0B 1.3X SMRTbell bead cleanup minimum required reagent.

### Protocol 3 of 10

#### 3.0A cDNA amplification

##### REAGENTS



Figure 3. 3.0A cDNA amplification minimum required reagent volume.

## Protocol 4 of 10

### 3.0B 0.9X cDNA SMRTbell bead cleanup

#### REAGENTS




Figure 4. 3.0B 0.9X cDNA SMRTbell bead cleanup minimum required reagent volumes.

## Protocol 5 of 10

### 4.0A Kinnex PCR

#### REAGENTS




Figure 5. 4.0A Kinnex PCR minimum required reagent volumes.

## Protocol 6 of 10

### 4.0B Pooling & 1.05X SMRTbell bead cleanup

#### REAGENTS




Figure 6. 4.0B Pooling & 1.05X SMRTbell bead cleanup minimum required reagent volumes.

## Protocol 7 of 10

### 5.0A Kinnex array formation

#### REAGENTS




Figure 7. 5.0A Kinnex array formation minimum required reagent volumes.

## Protocol 8 of 10

### 5.0B 1X SMRTbell bead cleanup

#### REAGENTS



Figure 8. 5.0B 1X SMRTbell bead cleanup minimum required reagent volumes.

## Protocol 9 of 10

### 6.0A Nuclease treatment

#### REAGENTS




Figure 9. 6.0A Nuclease treatment minimum required reagent volumes.

## Protocol 10 of 10

### 6.0B Final SMRTbell cleanup

#### REAGENTS



Figure 10. 6.0B Final SMRTbell cleanup minimum required reagent volumes.

# Consumables

| Supplier                 | Part Name  | Part Number   | Number Required | Note  |
|--------------------------|--|---------------|-----------------|---|
| SPT Labtech              | 40mm Upper Deck Riser                                      | 3276-01838    | 6               |   |
| SPT Labtech              | dragonfly® discovery Sterile LDV Reservoirs                | 4150-07203    | 15              | Reservoir types needed are dependent on the number of columns processed |
| SPT Labtech              | dragonfly® discovery Sterile Reservoirs                    | 4150-07204    | 23              |   |
| SPT Labtech              | dragonfly® discovery Sterile Syringes                      | 4150-07201    | 20              |   |
| SPT Labtech              | dragonfly® discovery Sterile, Ultra Low Retention Syringes | 4150-07209    | 18              | Number required depends on the number of columns processed              |
| SPT Labtech              | firefly® Strip Tips, 100µl, with Filters, Sterile          | 125-008-FF-FS | 25              |   |
| SPT Labtech              | firefly® Strip Tips, 35µl, with Filters, Sterile           | 050-008-FF-FS | 8               | Waste plate   |
| Alpaqua Engineering      | Alpaqua Magnum FLX   | A000400       | 1               |   |
| Thermo Fisher Scientific | Fisherbrand 1ml Deep Well                                  | 236600        | 4               |   |
| Bio-Rad                  | Hard Shell Plate (HSP)                                     | HSP-9601      | 1               |   |
| Eppendorf                | twin.tec PCR   | 30128648      | 4               |   |

**Table 3.** Consumables & labware required for Kinnex full-length RNA kit Library Prep on firefly.

## Protocol overview

This method was developed with an EZ-load 6 head genomics (v1.9.1 software) firefly using 24 samples, skirted twin.tec plates (30128648) and the Alpaqua Magnum FLX magnet. Use of alternative firefly configurations or labware may require further optimization.

This suite of protocols has been published with sample numbers assuming a 24 RNA sample input, with 4-sample multiplexing (3.3 Pooling Barcoded cDNA), resulting in 6 Kinnex PCR Pools & 48 Kinnex PCR samples (4.0 Kinnex PCR). Sample inputs may be adjusted using “Number of Samples” variables embedded in each protocol.

LDV reservoir dead volumes have been updated to 55 µL to account for limiting reagents where possible.


## Protocol 1 of 10 2.0A cDNA Synthesis

This protocol performs sections 2.1 – 2.3 of PacBio’s Kinnex Library Prep Kit (103-238-700 REV06).

### Prior to executing this protocol:

- **1.1** Measure RIN with Agilent 2100 Bioanalyzer Instrument using the RNA 6000 Nano kit
- **2.1.1** Thaw cDNA synthesis components as directed in protocol document

This protocol is compatible with 8 - 48 samples as written and has been published in the v1.8.6 firefly software.



**Figure 11.** 2.0A cDNA Synthesis deck layout.

## Protocol 2 of 10

### 2.0B 1.3X SMRTbell bead cleanup

This protocol performs section 2.4 of PacBio's Kinnex full-length RNA Kit (103-238-700 REV06).

#### Prior to executing this protocol:

- Equilibrate SMRTbell cleanup beads at room temperature for 30 minutes prior to use

This protocol is compatible with 8 - 48 samples as written and has been published in the v1.8.6 firefly software. To process > 48 samples, update 80% EtOH reservoir asset definition and EtOH aspiration steps to use four syringes.




Figure 12. 2.0B 1.3X SMRTbell bead cleanup deck layout.

## Protocol 3 of 10

### 3.0A cDNA amplification

This protocol performs section 3.1 of PacBio's Kinnex full-length RNA Kit (103-238-700 REV06).

#### Prior to executing this protocol:

- Aliquot Iso-Seq Primers to PCR Plate
- **3.1.1** Prepare reaction mix 3 (Iso-Seq cDNA PCR Master Mix)

This protocol is compatible with 8 - 96 samples as written and has been published in the v1.8.6 firefly software. To process > 48 samples, update 80% EtOH reservoir asset definition and EtOH aspiration steps to use four syringes.




Figure 13. 3.0A cDNA amplification deck layout.

## Protocol 4 of 10

### 3.0B 0.9X cDNA SMRTbell bead cleanup

This protocol performs section 3.2 of PacBio's Kinnex full-length RNA Kit (103-238-700 REV06).

#### Prior to executing this protocol:

- Equilibrate SMRTbell cleanup beads at room temperature for 30 minutes prior to use

This protocol is compatible with 8 - 96 samples as written and has been published in the v1.8.6 firefly software. To process > 48 samples, update 80% EtOH reservoir asset definition and EtOH aspiration steps to use four syringes.




Figure 14. 3.0B 0.9X cDNA SMRTbell bead cleanup deck layout.

## Protocol 5 of 10

### 4.0A Kinnex PCR

This protocol performs section 4.1 of PacBio's Kinnex full-length RNA Kit (103-238-700 REV06).

#### Prior to executing this protocol:

- **3.1.1** Aliquot Iso-Seq Primers to PCR Plate. Iso-Seq primer bc01-12 will be added to each sample individually and should not be added to the master mix
- **3.1.1** Prepare reaction mix 3 (Iso-Seq cDNA PCR Master Mix)

This protocol is compatible with 8 - 96 samples as written and has been published in the v1.8.6 firefly software. To process > 48 samples, update 80% EtOH reservoir asset definition and EtOH aspiration steps to use four syringes.




Figure 15. 4.0A Kinnex PCR deck layout.

## Protocol 6 of 10

### 4.0B Pooling & 1.05X SMRTbell bead cleanup

This protocol performs section 4.2 of PacBio's Kinnex full length RNA Kit (103-238-700 REV06).

#### Prior to executing this protocol:

- **4.2.1** Add exactly 23  $\mu$ L from each of the 8 PCR reactions into a 1.5 mL tube for a total volume of 184  $\mu$ L. An equal volume of each PCR product is necessary for efficient array assembly. Equilibrate SMRTbell cleanups beads at room temperature for 30 minutes prior to use
- Protocol deviation: Split each pool into 2 pools equaling 92  $\mu$ L for cleanup

This protocol is compatible with 8 - 96 samples as written and has been published in the v1.8.6 firefly software. To process > 48 samples, update 80% EtOH reservoir asset definition and EtOH aspiration steps to use four syringes.




Figure 16. 4.0B Pooling & 1.05X SMRTbell bead cleanup deck layout.

## Protocol 7 of 10

### 5.0A Kinnex array formation

This protocol performs section 5.1 of PacBio's Kinnex Library Prep Kit (103-238-700 REV06).

#### Prior to executing this protocol:

- **5.1.0** Place the DNA damage repair mix on ice at all times, and immediately return the DNA damage repair mix back to the freezer (-20°C) after use. Improper storage and handling of the DNA damage repair mix may result in poor library performance and should not be used for subsequent reactions
- **5.2.1** Prepare RM Master Mix(s)

This protocol is compatible with 8 - 96 samples as written and has been published in the v1.8.6 firefly software. To process > 48 samples, update 80% EtOH reservoir asset definition and EtOH aspiration steps to use four syringes.




Figure 17. 5.0A Kinnex array formation deck layout.

## Protocol 8 of 10

### 5.0B 1X SMRTbell bead cleanup

This protocol performs section 5.2 of PacBio's Kinnex full-length RNA Kit (103-238-700 REV06).

#### Prior to executing this protocol:

- Equilibrate SMRTbell cleanups beads at room temperature for 30 minutes prior to use

This protocol is compatible with 8 - 96 samples as written and has been published in the v1.8.6 firefly software. To process > 48 samples, update 80% EtOH reservoir asset definition and EtOH aspiration steps to use four syringes.




Figure 18. 5.0B 1X SMRTbell bead cleanup deck layout.

## Protocol 9 of 10

### 6.0A Nuclease treatment

This protocol performs section 6.2 of PacBio's Kinnex full-length RNA Kit (103-238-700 REV06).

#### Prior to executing this protocol:

- 6.1.1 Prepare Nuclease Master Mix

This protocol is compatible with 8 - 96 samples as written and has been published in the v1.8.6 firefly software. To process > 48 samples, update 80% EtOH reservoir asset definition and EtOH aspiration steps to use four syringes.




Figure 19. 6.0A Nuclease treatment deck layout.

## Protocol 10 of 10

### 6.0B Final SMRTbell cleanup

This protocol performs section 6.2 of PacBio's Kinnex full-length RNA Kit (103-238-700 REV06).

#### Prior to executing this protocol:

- 6.1.1 Prepare Nuclease Master Mix

This protocol is compatible with 8 - 96 samples as written and has been published in the v1.8.6 firefly software. To process > 48 samples, update 80% EtOH reservoir asset definition and EtOH aspiration steps to use four syringes.




Figure 20. 6.0B Final SMRTbell cleanup deck layout.